

iGEM TU/e 2014

Biomedical Engineering Eindhoven

University of Technology Room: Ceres 0.04 Den Dolech 2, 5612 AZ Eindhoven The Netherlands Tel. no. +31 50 247 55 59 2014.igem.org/Team:TU_Eindhoven

Date 14 October 2014

Antifouling



Table of contents

Title Reacting cells with DBCO-PEG10kDa	1	Stock solutions	3
	2	Protocol	3
	2	Dilution of Antibodics	2

1 Stock solutions

300 μM DBCO-PEG in DMSO: 250 μg + 83.3 μL DMSO

2.00 μM Anti-HA antibody TAMRA

Buffer: PBS with 0.1% BSA

2 Protocol

Prepare following tubes: (6 in total)

• Make 3 of each tube type.

Tube	[DBCO]	Cells (10 ⁷)	DBCO volume to add (μL) (300 μM)	DBCO/tag ratio
1	0	200 μL		
2	30 µM	200 µL	22.2	300.5

- Make sure you vortex the cell well before and after adding the DBCO-PEG
- React DBCO tubes for 1h to 6h in shaking block at 4°C and 300rpm
- Prepare FACS samples after 1h and 6h:
- Spin down the cells for 5 min at 13,400 rpm and discard the supernatant
- Resuspend with 1 mL ice cold PBS-0.1%BSA
- Spin down the cells for 5 min at 13,400 rpm and discard the supernatant
- Add 200 μL ice cold PBS-0.1%BSA

3 Dilution of Antibodies

Prepare following tubes: (Per DBCO 10kDa tube)

Tube	[Ab]	Antibody volume to add (μL) 2.00 μM	DBCO?
1	-		
2	-		Yes
3	316 nM	33.82	
4	100 nM	10.19	
5	316 nM	33.82	Yes
6	100 nM	10.19	Yes

- Let the reaction happen for 1 hour
- Prepare FACS samples:
 - o Spin down the cells for 5 min at 13,400 rpm
 - o Discard the supernatant and resuspend with 1 mL ice cold PBS-0.1%BSA
 - o Spin down the cells for 5 min at 13,400 rpm
 - o Discard the supernatant and put the pellets on ice until FACS
 - o Put on ice until FACS
 - O Right before FACS: resuspend with 200 μL ice cold PBS-0.1%BSA