Cloning procedure

Overview of a cloning procedure:

- 1. Digestion of a PCR product (= insert) and a vector
- 2. Dephosphorylation of the digested vector
- 3. Ligation of digested insert and vector
- 4. Transformation of E. coli
- 5. Screening of positive clones by control digestion and/or Colony-PCR
- 6. Linearization of plasmids for transformation of B. subtilis

1. Digestion using restriction enzymes

Calculation of a digestion:

$$\frac{\mathrm{U}\left(\mathrm{Enzyme}\right)}{\mu\mathrm{g}\left(\mathrm{x}\right)} = \frac{\mathrm{MW}\left(\mathrm{R}\right)}{\mathrm{MW}\left(\mathrm{x}\right)} \cdot \frac{\mathrm{n}\left(\mathrm{x}\right)}{\mathrm{n}\left(\mathrm{R}\right)}$$

with R Reference-DNA
x DNA to digest
MW Molecular weight
n Amount of cutting sites

- The phage Lambda (λ) serves as reference DNA with a MW of 48.5 kb
- The amount of λ 's cutting sites for a specific restriction enzyme is written down in the NEB catalog

EXAMPLE: digestion of 500 ng pXT with BamHI and HindIII

MW (pXT) = 6.892 kb
MW (
$$\lambda$$
) = 48.5 kb
n (pXT-BamHI) = 1
n (λ -BamHI) = 5
n (pXT-HindIII) = 1
n (λ -HindIII) = 6

BamHI:

→
$$\frac{\text{U (Enzyme)}}{\text{µg (pXT)}} = \frac{48.5 \text{ kb}}{6.892 \text{ kb}} \cdot \frac{1}{5} = 1.4$$

→
$$\frac{\text{U (Enzyme)}}{0.5 \,\mu\text{g (pXT)}} = \frac{1.4}{2} = 0.7$$

You need 0.7 U of BamHI to digest 500 ng of pXT

The concentration of the enzyme stock solution is 10 U/ μ l (for Fermentas enzymes!! \rightarrow NEB enzymes may have different units per μ l)

$$rac{10 U}{1 \mu l} = rac{0.7 U}{x}$$

$$\rightarrow$$
 x = 0.07 µl

To be really sure that almost everything of your DNA will be digested, double the calculated volume of your enzyme:

→
$$0.07 \, \mu l \cdot 2 = 0.14 \, \mu l$$

You need 0.14 µl of BamHI to digest 500 ng pXT.

HindIII:

$$ightharpoonup \frac{\text{U (Enzyme)}}{\text{µg (pXT)}} = \frac{48.5 \text{ kb}}{6.892 \text{ kb}} \cdot \frac{1}{6} = 1.2$$

⇒
$$\frac{\text{U (Enzyme)}}{0.5 \, \text{µg (pXT)}} = \frac{1.2}{2} = 0.6$$

You need 0.6 U of HindIII to digest 500 ng of pXT

The concentration of the enzyme stock solution is 10 U/μl

→
$$x = 0.06 \mu l$$

→
$$0.06 \mu l \cdot 2 = 0.12 \mu l$$

You need 0.12 µl of HindIII to digest 500 ng pXT.

Sample preparation:

- 500 ng pXT
- 0.5 μl *Bam*HI
- 0.5 μl HindIII
- X μl 10×Buffer
- X μl H₂O

Important: DO NOT exceed the volume of the enzymes above 1/10 of the final volume

→ Incubate for 1.5 to 2 h at 37°C (if you use FastDigest enzymes: incubate for 15 min at 37°C)

2. Dephosphorylation of digested vectors

Important: DO NOT dephosphorylate your inserts!

FastAP (Thermosensitive Alkaline Phosphatase; Fermentas):

- Add 1 μl of Fast AP to the digested vector
- Incubate at 37°C for 15 min

SAP (Shrimp Alkaline Phosphatase; Fermentas):

- Add 1 μl of SAP to the digested vector
- Incubate at 37°C for 30 min
- Add another 1 μl of SAP and incubate further 30 min at 37°C

After dephosphorylation: Purify the digestion samples by using the PCR Purification Kit

3. Ligation

Calculation of a ligation:

$$\frac{\text{ng (Vector)}}{\text{kb (Vector)}} = \frac{ng (Insert)}{kb (Insert)}$$

EXAMPLE: ligation of 200 ng pXT with a 1000 bp DNA fragment

$$\rightarrow$$
 x = 29 ng

Molecular ratio of vector (200 ng pXT) to insert (29 ng of a 1000 bp DNA fragment) = 1:1

In the most cases you will need a 3-times molar excess of insert to vector

→ 29 ng · 3 =
$$87 \text{ ng}$$

Sample preparation:

- 200 ng pXT
- 87 ng Insert
- 1 µl T4 DNA Ligase (Fermentas)
- X μl 10×Ligase-Buffer (Fermentas)
- X μl dH₂O

- → Keep the final volume as small as possible!
- → Incubate for 1-3 h at room temperature

4. Transformation of competent E. coli cells (CaCl₂)

- Thaw 100 µl of competent cells on ice
- Add ½ volume of your ligation sample to the cells
- Incubate for 20 min on ice
- Heat shock the cells at 42°C for 1.5 min
- · Cool down the cells for 2 min on ice
- Add 400 µl LB to the cells
- Shake the cells for 1 h at 37°C
- Plate 50 μl of the cell solution on a LB+antibiotic agar plate
- Centrifuge the remaining cell solution (1 min, 13 krpm, room temperature)
- Discard the medium and resuspend the pellet in the remaining medium drop
- Plate the entire resuspended cells on a LB+antibiotic agar plate
- Incubate the plates at 37°C for overnight

6.1 Control digestion

Sample preparation:

- 200 ng plasmid
- 0.3 μl enzyme 1
- 0.3 μl enzyme 2
- Appropriate buffer
- dH₂O
- → Incubate for at least 1 h at the dedicated temperature
- → Mix with DNA loading dye and load the entire sample on an agarose gel

6.2 Colony-PCR → See PCR methods

7. Linearization of integrative plasmids for transformation of B. subtilis

Sample preparation: (30 µl final volume)

- 25 μl plasmid (from Mini Prep)
- 1 µl enzyme
- 3 µl buffer
- 1 μl dH₂O
- → Incubate for 1.5-2 h at 37°C
- → Use 15 µl of the sample for *B. subtilis* transformation

Examples for linearization enzymes of different vectors: (Check in advance that the indicated enzymes do not cut your insert!)

Vector	Enzyme for linearization
рХТ	Scal
pAC5	PstI
pAC6	PstI or Scal
pAC7	PstI
pSWEET	Scal
pALFLAGrsiW-1	Pvul
pDG1662	Scal
pDG1663	Ncol
pDG1731	Scal

Protocol generously provided by the lab Prof. Thorsten Mascher Großhadernerstr. 2-4 82152 Planegg-Martinsried www.syntheticmicrobe.bio.lmu.de