



University of Melbourne iGEM 2014 Lab Procedure

Procedure	Name:	PCR Protocol		
	Description:	For amplification of DNA. Taken from the Q5 High-Fidelity 2X Master Mix instructions. <a href="#">PCR Video</a>		
	Trigger:			
Last updated	Name:	Elizabeth Brookes	Date:	15.08.14
You will need	Time:	2 hours		
	PPE:	Standard equipment (lab coat, gloves, etc.)		
	Equipment:	<ul style="list-style-type: none"> <li>○ Thermocycler/PCR machine</li> <li>○ Ice box</li> </ul>		
	Materials:	<ul style="list-style-type: none"> <li>○ 4 x PCR tubes</li> <li>○ 1mL of Q5 High-Fidelity 2X Master Mix</li> <li>○ 10μL of 10μM Forward Primer</li> <li>○ 10μL of 10μM Reverse Primer</li> <li>○ 8μL of Template DNA</li> <li>○ 72μL of Nuclease-Free Water</li> </ul>		
Step 1	<p>For a standard 50μL reaction, mix the following into a PCR tube <b>on ice</b> and repeat for all 4 tubes:</p> <ul style="list-style-type: none"> <li>○ 25μL of Q5 High-Fidelity 2X Master Mix</li> <li>○ 2.5μL of 10μM Forward Primer</li> <li>○ 2.5μL of 10μM Reverse Primer</li> <li>○ 2μL of Template DNA</li> <li>○ 18μL of Nuclease-Free Water to make up to 50μL</li> </ul>			
Step 2	Gently mix the reagents. If required, collect all liquid to the bottom of the PCR tubes with a quick spin.			
Step 3	Place the PCR tubes in the PCR machine and close the lid.			
Step 4	<p>Program the PCR machine:</p> <ol style="list-style-type: none"> <li>1. Press &lt;Create&gt;.</li> <li>2. Press &lt;Start&gt;.</li> <li>3. Using the arrow keys and number pad, program the following temperatures and times: <ul style="list-style-type: none"> <li>• Initial denaturation: 98°C, 5 minutes.</li> <li>• 30 cycles of: <ul style="list-style-type: none"> <li>• Denaturation: 98°C, 30 seconds.</li> <li>• Annealing: 57°C, 30 seconds.</li> <li>• Extension: 72°C, 30 seconds.</li> </ul> </li> <li>• Final extension: 72°C, 5 minutes.</li> <li>• Hold: 4°C, infinite time.</li> </ul> </li> <li>4. Press &lt;Start&gt;.</li> <li>5. Enter reaction volume (50μL) when prompted.</li> <li>6. Press &lt;Start&gt; to begin the reaction.</li> </ol>			
Step 5	Wait for the PCR machine to finish (approximately 40 minutes). The screen will show which stage the machine is up to.			
Step 6	Collect the tubes from the machine. If they are not to be used immediately, place in the -20°C freezer. Otherwise, place in ice.			

Version history	<i>V1 replaced by the Q5 High-Fidelity 2X Master Mix protocol to create V2.</i>
-----------------	---