

Golden Gate Cloning

- Add 100 ng of the linearized vector backbone and equimolar amounts of the other assembly pieces to 15 µl total volume assembly reaction mixture as follows

Linearized vector backbone (100ng)

+		each additional assembly piece (to equimolar with backbone)
+	1,5 µl	10X NEB T4 Puffer
+	0,15 µl	100XBSA*
+	1 µl	Bsal
+	1 µl	NEB T4 Ligase, 2 million cohesive end units/mL
+		dH ₂ O to

15 µl

NOTE: It is essential to use a High Concentration Ligase

* Bsal is only 10 % active at 37 °C without the addition of BSA.

- Perform the assembly reaction in a thermocycler as follows: (Engler 2009)

3	min	@ 37 °C }
4	min	@ 16 °C } 25 cycles
5	min	@ 50 °C }
5	min	@ 80 °C } 1 cycle

- Transform into competent *E.coli*

Protocol generously provided by the lab
Prof. Thorsten Mascher
Großhadernerstr. 2-4
82152 Planegg-Martinsried
www.syntheticmicrobe.bio.lmu.de