



University of Melbourne iGEM 2014 Lab Procedure

Procedure	Name:	Restriction Digestion		
	Version:	3		
	Description:	Restriction digestion allows us to excise our gene of interest from a plasmid and linearise vector plasmids for ligation. Restriction Digestion Video		
Last updated	Name:	Robyn Esterbauer	Date:	22.06.14
You will need	Time:	1 ½ hours		
	PPE:	Lab coat Safety Glasses Gloves		
	Equipment:	P20 Pipette and tips Sample tubes and rack Ice and container		
	Materials:	Restriction digestion enzyme (e.g. Not1) Restriction digestion enzyme buffer (e.g. CutSmart Buffer) Vector Milli-Q Water Bovine Serum Albumin (BSA)		
Step 1	Aliquot 30µL of each vector being digested (i.e. pMK-RGD, pMK-Mag1, pET21d, pGEX6P-3) into separate microfuge tubes. Label the tubes appropriately eg. RD-pET-Mag1 22.06.14 RE iGEM			
Step 2	<p>Pipette into each sample tube, to make a total 50µL (note, you should actually follow the instructions on the NEB restriction enzyme cards):</p> <ul style="list-style-type: none"> ○ 12µL Milli-Q water ○ 5µL restriction digestion enzyme buffer (e.g. CutSmart Buffer) ○ 1-2µL restriction digestion enzyme (e.g Not1) <p><i>Note: Adding the Milli-Q water first means that you can use the same pipette tip for each aliquot. For all following steps use a different tip to avoid cross-contamination.</i></p> <p><i>Note: Always add the restriction digestion enzyme last and leave it in the freezer until you are ready to use it.</i></p> <p>Note that you can decrease the total reaction volume if you want the concentration of DNA to be higher. Note that there is a spreadsheet on the dropbox in the iGEM protocols folder (search for “Restriction digest calculator.xlsx”) which makes it easy to calculate new reaction volumes</p> <p>For double digestions, the concentration of restriction enzymes should not be greater than 10% by volume of the total reaction volume (see the NEB Double Digestion site, subheading “Setting up a Double Digestion”). For example, if the total reaction volume is 50 µL, then don’t add more than 5 µL of restriction enzymes.</p> <p>Also, note that we have limited amounts of SpeI; so, consider adding less SpeI and incubating for longer.</p>			
Step 3	Mix the contents in the tubes using pipette tips to move it up and down. Ensure that a different tip is used between tubes.			

Step 4	Spin using the desktop centrifuge on pulse for about 30 seconds to make sure all the contents are mixed with no bubbles
Step 5	<p>Place in 37°C incubator for at least 1 hour (don't follow the instructions on the NEB cards for the reaction time).</p> <p>Note that if you are preparing the sample for a ligation reaction, you really want a thorough digest and cannot accept partial digestion of a parent plasmid. Thus, according to George, you should incubate for at least two hours. If partial digestion is acceptable, then you get away with much less time.</p>
Version history	<p>Version 3 created on 2014-07-30 by Sean Lowe.</p> <p>Version 2 created on 16.07.14 by Elizabeth Brookes</p> <ul style="list-style-type: none">○ More detail added to all sections.