

DNA Kit Plate Instructions

To use the DNA in the Distribution Kit you may follow these instructions:

1. With a pipette tip, punch a hole through the foil cover into the corresponding well of the Biobrick™-standard part that you want. [Make sure you have properly oriented the plate](#). We recommend that you do not remove the foil cover, as it could lead to cross contamination between the wells.
2. Pipette 10uL of dH₂O (distilled water) into the well. Pipette up and down a few times and let sit for 5 minutes to make sure the dried DNA is fully resuspended. We recommend that you do not use TE to resuspend the dried DNA.
3. [Transform](#) 1 or 2uL of the resuspended DNA into your desired competent cells, plate your transformation with the appropriate antibiotic* and grow overnight.
4. Pick a single colony and inoculate broth (again, with the correct antibiotic) and grow for 16 hours.
5. Use the resulting culture to [miniprep](#) the DNA AND make your own glycerol stock (for further instruction on making a glycerol see [this page](#)). We recommend using the minipreped DNA to run QC tests, such as restriction digests and sequencing.

** To know which antibiotics to use, look at the plasmid that the part is in. The [naming scheme](#) for plasmids is specifically designed to indicate antibiotic resistance.*

Note: There is not enough DNA in each well to perform anything but [transformations](#)