



J23101	307	8.2	5	1	-	1	-	34.8
E0240	81	31	5	-	1	-	1	12.0
pSB1A2	181	13.8	5	1	-	-	1	29.2

## 2nd Day – August 30th, 2014

### Agarose (1%) gel electrophoresis

- Total digestion product volume of each part DNA was applied to the agarose gel 1% stained with ethidium bromide (0.5 ug/mL) and bands showing the corresponding size were excised from the gel and purified according the Wizard SV gel and PCR clean up kit protocol (Promega).

- The purified products were quantified using Nanodrop™ 1000 (Thermoscientific)

### Ligation

#### Assembly 3

<b>Linear plasmid (50 ng)</b>	<b>2 uL</b>	
Insert: Plasmid 6:1	J23115	E0240
	2 uL	10 uL
T4 DNA Ligase buffer (10X)	2 uL	
T4 DNA ligase (1U)	1 uL	
H <sub>2</sub> O to 20 uL	5 uL	

#### Assembly 2

<b>Linear plasmid (50 ng)</b>	<b>2 uL</b>	
Insert: Plasmid 6:1	J23101	E0240
	2 uL	10 uL
T4 DNA Ligase buffer (10X)	2 uL	
T4 DNA ligase (1U)	1 uL	
H <sub>2</sub> O to 20 uL	5 uL	

To determinate the amount of DNA necessary we used the following equation

$$\text{Insert ng} = \text{plasmid ng} \times \frac{\text{insert bp}}{\text{plasmid bp}} \times \text{insert: plasmid ratio}$$

The Ligation reactions were Incubated overnight at 37 °C.

Meanwhile, we prepared and sterilized in the autoclave tubes of liquid LB medium and glycerol 40%.

### **3<sup>rd</sup> Day - August 31th, 2014**

Transformation in competent *Escherichia coli* TOP10 (Invitrogen) using cloranphenicol (70 ug/mL) for selection of positive clones with Assembly 1 and ampicillin (100 ug/mL) ) for selection of positive clones with Assembly 2. The Transformation protocol is presented below:

#### **Transformation in Escherichia coli (TOP10 - Invitrogen)**

##### **Materials**

- 1,5ml tube
- Styrofoam box with ice
- Agar plate with antibiotic
- Competent cells
- Plasmidial DNA
- Centrifuge
- Water bath at 42°C
- Liquid LB medium
- Shaker

##### **Method**

- Briefly spin the competent cells and put then on ice
- Add 50ng of plasmidial DNA in a 1,5ml tube
- Add the 50ul of competent cell in the same tube
- Keep the tube on ice for 25min
- Put the tube in a 42° water bath for 2 min
- Put the tube on ice for 5 min
- Add 200ul of liquid LB
- Incubate at 250rpm/37°C/1 hour
- Plate the the solution in a Agar plate with the appropriate antibiotic
- Incubate the plate at 37°C overnight

### **4<sup>th</sup> Day - September 1st, 2014**

Inoculate 2 colonies in 5 mL LB mediumcontaining the correspondent antibiotic.

### **5<sup>th</sup> Day - September 2nd, 2014**

We prepared glycerol stocks of the clones, using 250 uL of clone culture + 250 uL of sterile glycerol 40%. The samples were stocked at -80°C.

4,75 mL of each clone culture was submitted to Plasmidial Extraction using PureLink Invitrogen Protocol.


### Assemble confirmation - September 3rd, 2014

We performed an EP Digestion as follows to confirm the assemble formation:

	Volume to 300 ng (uL)	Buffer x10 (uL)	FastDigest® EcoRI NEB (uL)	FastDigest® PstI NEB (uL)	H <sub>2</sub> O to 10 uL (uL)
<b>Assembly 3</b>	3	1	0.5	0.5	5
<b>Assembly 2</b>	3	1	0.5	0.5	5

Agarose (1%) gel electrophoresis

- Confirmation reaction product for one clone was applied to the agarose gel 1% stained with ethidium bromide (0.5 ug/mL). The result is shown in the figure bellow:



The results confirm Assembly 3 and Assembly 2 formation.