



University of Melbourne iGEM 2014 Lab Procedure

| | | | | |
|---|---|--|-------|----------|
| Procedure | Name: | Mini-purification of His-tagged proteins | | |
| | Version: | 2 | | |
| | Description: | How to purify His-tagged proteins using batch purification. | | |
| | Trigger: | | | |
| Last updated | Name: | Sean Lowe | Date: | 29.07.14 |
| You will need | Time: | | | |
| | PPE: | Gloves Lab coat | | |
| | Equipment: | Liquid nitrogen Dewar Plastic double walled well Ultracentrifuge pH Meter Peristaltic Pump Pipette and tips | | |
| | Materials: | Frozen cell pellets Ni-NTA beads Lysis Buffer 100mL pH 8 <ul style="list-style-type: none"> ○ 20mM Tris ○ 2mM Benzamidine ○ 0.1mg/mL Soybean Trypsin Inhibitor ○ 10% Glycerol ○ 1% Triton x100 (if required) ○ 300mM NaCl ○ 5mM Imidazole Lysozyme Iodoacetamide DNase Equilibration buffer Bradford reagent Base Buffer 500mL pH 8 <ul style="list-style-type: none"> ○ 20mM Tris ○ 2mM Benzamidine ○ 10% Glycerol | | |
| A) Prepare buffers | | | | |
| Step 1 | Prepare lysis buffer according to recipe listed above. | | | |
| Step 2 | Add lysozyme to the lysis buffer to a final concentration of 0.1 mg/mL. Add iodoacetamide now, if required. | | | |
| B) On The Day: Lysis using freeze-thaw in Eppendorf tubes | | | | |
| Step 1 | Scrape off a bit of the frozen cell pellet into the lysis Eppendorf tube | | | |
| Step 2 | Add the prepared lysis buffer to the cell pellet in a ratio of 1ml:1g wet cell pellet. | | | |
| Step 3 | Freeze-thaw by dunking the lysis Eppendorf tube into liquid nitrogen, then into a 37 degree water bath. Repeat another 2 times. | | | |

| | |
|---|--|
| Step 4 | Add DNase to the lysis Eppendorf tube at a concentration of 10 ug/mL. Add 5 mM MgCl ₂ (3.75 uL of a 2 M stock for 1.5 mL of lysis buffer) |
| Step 5 | Incubate with the DNase for 15 min on the rotary suspension mixer at room temperature |
| Step 6 | Centrifuge spin at 4 C at 12,000 rpm for 30 minutes |
| Step 7 | Aliquot out 40µL of sample for running on a gel |
| C) While tubes are spinning in step B6: Preparation the Ni-NTA beads: | |
| Step 1 | Vortex the slurry that the beads come in |
| Step 2 | Add an appropriate amount of Ni-NTA resin to an eppendorf tube. (minimum 100µL) |
| Step 3 | Centrifuge tube for 2 minutes at 700 × g and carefully remove and discard the supernatant. |
| Step 4 | Add two resin-bed volumes of Equilibration Buffer (250 uL) and mix until the resin is fully suspended. |
| Step 5 | Centrifuge tube for 2 minutes at 700 × g and carefully remove and discard Equilibration buffer. |
| D) Sample addition and incubation | |
| Step 1 | Use a pipette to carefully take up the supernatant from the lysis Eppendorf tube. Save the insoluble pellet remaining for SDS-PAGE analysis later. |
| Step 2 | Add the prepared protein extract to the tube from C5 and mix on an end-over-end rotator for at least 30 minutes in the cold room |
| E) Wash to remove nonspecific binding | |
| Step 1 | Centrifuge the tube for 2 minutes at 700 × g. Save supernatant for downstream analysis, labelling “Unbound” |
| Step 2 | Wash the resin with two resin-bed volumes (2x300µL) of Wash Buffer. Centrifuge the tube for 2 minutes at 700 × g. Save supernatant for downstream analysis. |
| Step 5 | Repeat step E2 using a new collection tube until there is no more protein in the washes as detected using the Bradford reagent |
| F) Elution | |
| Step 1 | Elute bound His-tagged proteins using one resin-bed volume of Elution Buffer (300 uL). Centrifuge tubes for 2 minutes at 700 × g. Carefully remove and save the supernatant. Repeat this step at least twice, saving each supernatant fraction in a separate tube. |
| Step 2 | Monitor protein elution by measuring the absorbance of the fractions at 280nm or by Coomassie Plus (Bradford) Assay Reagent |
| Step 3 | The eluted protein can be directly analyzed by SDS-PAGE. |
| Version history | <i>Describe the changes made in each new version of the protocol here.</i> |