

Todd-Hewitt-Broth (THB)

•

Beef heart infusion	3,1 g
Peptone	20 g
Glucose	2,0 g
NaCl	2,0 g
Na ₂ HPO ₄	0,4 g
Na ₂ CO ₃	2,5 g
H ₂ O	ad.1000 ml

➤ autoclave

Protocol generously provided by the lab
Prof. Thorsten Mascher
Großhadernerstr. 2-4
82152 Planegg-Martinsried
www.syntheticmicrobe.bio.lmu.de